

Microorganisms and Biofilms

LS.1.002

Electron and light microscopy of yeast biofilm

K. Dobranska¹, F. Ruzicka², J. Nebesarova³, Z. Burdikova⁴, J. Dluhos⁵, J. Collakova⁶, R. Chmelik⁶
O. Samek¹, V. Krzyzanek¹

¹ISI ASCR, Brno, Czech Republic

²St. Anne's hospital, Brno, Czech Republic

³Biology center of the ASCR, Ceske Budejovice, Czech Republic

⁴Institute of physiology of the ASCR, Praha, Czech Republic

⁵Tescan, Brno, Czech Republic

⁶UT Brno, Brno, Czech Republic

dobranska@isibrno.cz

Keywords: biofilm, cryo-SEM, light microscopy

Microbes like bacteria or yeast have different ways of living; they are able to live like planktonic cells or in more dangerous form in adherence to surfaces or interfaces, where are embedded in a matrix of extracellular polymeric substances that they can produce. Biofilm allows protection for the microbial cells from attacks by the immunity system as well as from the effect of antibiotics [1]. Therefore, study of biofilms is important for clinical research. It may help to develop more efficient treatment strategies for biofilm infections. Here we investigate the extracellular matrix of yeast *Candida parapsilosis* and *Candida albicans*.

Yeast and its biofilm have been investigated by scanning electron microscopy (SEM) as well as light microscopy. For surface imaging of the samples both classical and cryo-SEM techniques were employed and compared. For structural characterization also focused ion beam SEM (FIB-SEM) and the cryo-SEM freeze-fracturing technique were applied. FIB-SEM was used for both precise cross section preparations as well as for tomographical acquisition of 3D dataset using the signal of backscattered electrons (BSE) [2]. Freeze fracture provides planar views of the internal organization of membranes or biofilms and thus gives unique structural information.

For in vitro imaging coherence-controlled holographic microscopy (CCHM) [3] is employed; the main advantage remains the quantitative phase contrast imaging for non-invasive label-free live cell. Another light microscopy technique is the two-photon fluorescent confocal microscopy that gives much higher resolution than classical confocal light microscopes. Various labelling are used and compared.

A comprehensive view of the structure of fully hydrated system of yeast cultures living in the polymer matrix is obtained due to using a number of different imaging techniques. The sample preparation is the most critical point and therefore it is important and valuable to compare different microscopic approaches with different sample handling.

1. R.M. Donlan, and J.W. Costerton, Clin. Microbiol. Rev., 15, 167-193 (2002).
2. M. Zdražil in "The Step Towards an Ultimate Multifunctional Tool for Nanotechnology" et al., (3M Nano 2012 Proceedings, X'ien) (2012) 175
3. P. Kolman, and R. Chmelik, Optics Express, 18 (2010), p. 21990
4. The authors acknowledge the support by the grants CZ.1.07/2.3.00/20.0103 and CZ.1.05/1.1.00/02.0068 (EC and MEYS CR), P205/11/1687 (GACR), TE01020118 (TACR), FR-TI4/660 (MIT CR) and Z60220518 (ASCR).




Figure 1. (A) Surface of yeast *C. parapsilosis* grown on the agar (cryo-SEM); (B) Freeze-fracturing and cryo-SEM of grown yeast culture of *C. parapsilosis* on cover glass; (C) Maximal projection of *C. parapsilosis* grown on cover glass with Calcofluor labelling; (D) FIB-SEM tomograph of *C. albicans* after freeze-substitution; (E,F) Series of quantitative phase contrast (CCHM) of growing yeast *C. parapsilosis* culture at times 0h and 6h.